

HOLLYWOOD SPECTRA™

LUTRONIC
Intelligent Care

Official platform for the
Hollywood Laser Peel®

Highest peak power
for efficacy & virtually
no downtime

Proven aesthetic
performance & versatility


THE WORLD'S BEST JUST GOT BETTER

Dual wavelength **1064/532nm*** platform safely treats **all skin types and numerous aesthetic concerns****

HyperSurge™ Resonator delivers **high peak acoustic pressure** that exceeds most picosecond lasers and **dual pulse modes** for treatment versatility

- 5ns for effective, no downtime photomechanical treatment of unwanted pigments
- 0.3ms for controlled photothermal skin revitalization
- Uniquely enables the true Hollywood Laser Peel and much more

With **8000+ systems installed**, it has globally proven reliability, ease-of use, and unmatched versatility for practice ROI.



PROPRIETARY INTELLIBEAM™
radial micro-beam homogenizer
ensures energy consistency



**1st US FDA
Cleared
Device to Treat
Melasma**

HIGH PERFORMANCE, QUALITY & ROI FOR YOUR PRACTICE

Hollywood Laser Peel®

Exfoliates Skin

Reduces Pigmentation & Sun Damage

Helps Stimulate Collagen

Softens Skin Tone & Texture

Fast Treatments

No Discomfort or Downtime



Check the buzz: [#hollywoodlaserpeel](https://twitter.com/hollywoodlaserpeel)

Better Manage Melasma

Successfully treat unwanted pigmentation and stimulate collagen in all skin types.

BEFORE & AFTER 7 TREATMENTS



RESULTS MAY VARY. COURTESY OF MELANIE PALM, MD

“The Hollywood Spectra platform provides leadership performance in reducing unwanted pigment and stimulating collagen with best in class laser toning. With peak power comparable to pico systems but with much lower operating costs, Spectra is a reliable and versatile workhorse laser that we use every day. And the patented Hollywood Laser Peel is a fast and popular no downtime treatment for many patients seeking the red-carpet look.”

OMAR IBRAHIMI, MD, PHD

*Optional 585/650nm Dye Handpieces available to treat the less prevalent green & light blue tattoos.

**Laser toning, melasma, pigmented lesions, epidermal nevi, tattoo removal, acne scars, fine wrinkles, non-ablative resurfacing, skin revitalization, PIH, nevus of ota, lentigo, freckles, SKs.

19 Fortune Drive, Billerica, MA 01821 USA 888-588-7644 usa.lutronic.com contact@lutronic-usa.com

©2020. LUTRONIC Group of Companies. All rights reserved. Lutronic, its logo, Lutronic Intelligent Care, Hollywood Spectra, Hollywood Laser Peel, HyperSurge, and Intellibeam are trademarks or registered trademarks of the Lutronic Group of Companies. This product or the use of this product is covered by one or more US and foreign patents or patent applications pending. MARK_HSP_SLS_000

LUTRONIC
Intelligent Care™